



EVALUATION OF THE STERILITY OF SINGLE-USE VIALS UNDERGOING MULTIPLE ACCESS FOLLOWING APPLICATION OF A CLOSED SYSTEM TRANSFER DEVICE (CSTD)

Perks W, Carating H, Iazzetta J, Charbonneau LF, Walker SE.

Department of Pharmacy, Sunnybrook Health Sciences Centre and Leslie Dan Faculty of Pharmacy, University of Toronto.

BACKGROUND

- Closed system transfer devices (CSTD) such as Cyto-Set® (1), PhaSeal® (2), and Equashield® (3) are designed to reduce hazardous drug exposure from preparation to administration (4)
- While these devices protect staff and have been provided an ONB designation by the FDA, as a closed system they can also minimize microbial contamination because the devices are airtight and leak-proof, potentially preventing microbial ingress.
- Equashield® was rated favourably in terms of ease of use and therefore used as the chosen CSTD for this study (4,6)

OBJECTIVES


We undertook this study to test whether attaching Equashield® adaptors to simulated single-use vials (Tryptic Soy Broth 20 mL vials):

- Could prevent or minimize microbial contamination and extend the “use-by” date following multiple withdrawals under extreme-use-conditions
- Could verify the results of previously published data (5) applied to real-world conditions with multiple staff members and multiple areas at Sunnybrook Health Sciences Centre Pharmacy

METHODS

- Each lot of TSB growth medium was secondarily tested as a positive control by inoculation with less than 10⁻² of *S. epidermidis* ATC 12228, incubation at 37°C and review at 24 and 48 hours afterwards
- As a negative control, an unopened vial of TSB 20 mL from each lot was incubated for the duration of the study

Weekly: 18 vials of 20mL Tryptic Soy Broth + Equashield® Adaptor



Each vial group was placed in one of two settings for 5 day period:
1. BSC (biological safety cabinet Class II Type B2) running 24/7 (3 A/B/C vials)
2. Laminar airflow hood turned off nightly (3 A/B/C vials)

- At the end of day 5 each week, the vials were collected, incubated at 37°C for 14 days and inspected visually every 2 days for microbial growth
- This process was repeated for 16 weeks for a total of 288 vials (32 vials per group)

FIGURE 1: Study Design Procedure

DOWNLOAD THIS POSTER

Interested in a copy of this poster and other Sunnybrook Posters?

Go to: http://metrodis.org/SB_PPC2015.html

and download the poster from this site or scan the QR code →



RESULTS

Table 1: Distribution of TSB Vial Access and Positive Results (16 week period)

Area	AREA 1 (TWO BSCs RUNNING 24/7)			AREA 2 (ONE BSC RUNNING 24/7)			AREA 3 (THREE LAMINAR AIR FLOW HOODS, RUNNING DURING DAYTIME ONLY)		
	Vial A	Vial B	Vial C	Vial A	Vial B	Vial C	Vial A	Vial B	Vial C
Positive Results over number of accesses	0 / 32	0 / 160	0 / 320	0 / 16	0 / 80	0 / 160	0 / 48	0 / 240	0 / 480

- All positive control vials demonstrated growth within 48 hours
- All negative control vials showed no growth throughout the study
- All accessed vials remained sterile following storage at room temperature for 5 days and subsequent incubation for 14 days
- None of the 192 vials (B & C) accessed 1440 times showed evidence of contamination. The 95% confidence interval of the contamination rate is 0.000 to 0.035%
- None of 96 vials with CSTD attached but without broth removed demonstrated contamination

COST SAVINGS IMPLICATIONS

- Annual drug wastage in the cancer centre pharmacy, adhering to chemical stability and preparation in an ISO Class 5 environment, but not discarding partial vials after 6 hr as per USP<797> was ~\$185,000 in 2014/15. This represents ~1% of drug expenditures.
- Annual estimates of drug wastage in 2015/16, abiding USP<797> would be in excess of \$2.7M.
- If we were to use a CSTD on every single use vial in our system, the drug cost savings to our system is estimated to exceed \$2.5M annually.
- The incremental cost of the CSTD within the outpatient oncology program is estimated at \$400,000 annually, based on drug expenditures of \$20.4M and approximately 25,000 patient treatment visits.
- Net cost savings would be ~\$2.1M




FIGURE 2: (A) Negative Control Vial (B) Study vial with Equashield® adaptor attached

(C) Equashield® syringe

CONCLUSIONS

- Attachment of a CSTD adaptor to single-use vials within an ISO-5 environment has the ability to maintain sterility following multiple withdrawals
- These results were consistent when the vials were exposed to continuous ISO 5 air quality and when exposed to poorer than ISO 5 air quality (Laminar air flow hood off)
- These results remained consistent in the setting of multiple operators following training on the use of the CSTD

REFERENCES

1. Gebel J. Evaluation of the microbial tightness of closed system transfer devices by simulating airborne and touch contamination. *Critical Care*. 2015;19(Suppl 1):P85.
2. Edwards M, et al. Cost savings realized by use of the PhaSeal(R) closed-system transfer device for preparation of antineoplastic agents. *Journal of Oncology Pharmacy Practice*. 2013;19(4):338-347.
3. <http://www.capho.org/sites/default/files/nops/Evaluation%20of%20a%20New%20Closed%20System%20Transfer%20Device%20-%20Flay%20Charbonneau.pdf>
4. Vyas N, et al. Occupational exposure to anti-cancer drugs: A review of effects of new technology. *Journal of Oncology Pharmacy Practice*. 2013;20(4):278-287.
5. Rowe E, Savage S, Rutala W, Weber D, Gergen-Teague M, Eckel S. Economic and Microbiologic Evaluation of Single-Dose Vial Extension for Hazardous Drugs. *Journal of Oncology Pharmacy Practice*. 2012;8(4):e45-e49.
6. Heather Scott, Susan Singh, Flay Charbonneau. Evaluation of a New Closed System Transfer Device. *Canadian Association of Pharmacy in Oncology Conference*. St. John's, Nfld. May 21-24, 2015.